

CPNs™	Ex/Em maxima (nm)	Fluorescent colour	Spectrally comparable* fluorophores	Commonly used filter sets (Ex/Em)
CPN™ 1130	750 / 1130	NIR-II		
CPN™ 900	650 / 900	NIR-I		
CPN™ 680	400 / 680	Red		
CPN™ 610	480 / 610	Orange		
CPN™ 550	470 / 550	Yellow	FITC, Alexa Fluor® 488, GFP, YFP	475/70 BP 530/86 BP
CPN™ 510	400 / 510	Green B		
CPN™ 510	450 / 510	Green	FITC, Alexa Fluor®488, Dylight 488,GFP, YFP	475/70 BP 530/86 BP
CPN™ 475	390 / 475	Blue	AMCA, eBFP, DAPI, Hoechst 33342, Hoechst 33258	377/60 BP 447/60 BP
CPN™ 435	390 / 435	Indigo	AMCA, eBFP, DAPI, Hoechst 33342, Hoechst 33258, Alexa Fluor® 405	377/60 BP 447/60 BP
CPN™ 420	390 / 420	Violet	AMCA, eBFP, DAPI, Hoechst 33342, Hoechst 33258, Alexa Fluor® 405	377/60 BP 447/60 BP



LINKBRIGHT™

CPN: Thiol Oligonucleotide Conjugation Kit

Manual and protocol

60min reaction time, 7 easy steps

This product is for research use only and is not intended for diagnostic use.

LINKBRIGHT™ Conjugation Kits

are available for IgG antibody, oligonucleotide and protein linkage via Amine & Thiol in the above CPN wavelengths.

Purchaser Notification

These high-quality reagents and materials must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Read the Safety Data Sheet for each product available, other regulatory considerations may apply.

Obtaining Support

Search FAQ's at www.streambio.co.uk/FAQs or submit a question directly to Technical Support techsupport@streambio.co.uk

SDS

Safety Data Sheets (SDSs) are available at www.streambio.co.uk/resources-downloads/

Certificate of Analysis

The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on request

For Research Use Only. Not for use in diagnostic procedures.

STREAM BIO LTD AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL STREAM BIO LTD AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

Limited Warranty

Stream Bio Ltd and/or its affiliate(s) warrant their products as set forth in the Stream Bio' General Terms and Conditions of Sale found on Stream Bio's' website at www.streambio.co.uk/termsandconditions If you have any questions, please contact Stream Bio Ltd at info@streambio.co.uk or techsupport@streambio.co.uk

Materials supplied and storage

Vial / cap colour	Item (colour coded)	Quantity			Storage
1. Grey*	(150µl) LINKBRIGHT™ Thiol	1 vial	3 vials	10 vials	-20°C
2. Pink	Solution HT (200µl)	1 vial	1 vial	1 vial	-20°C
3. Orange	Solution SC (200µl)	1 vial	1 vial	1 vial	-20°C

When stored as directed, the kit components are stable for at least 6 months.
1x vial of CPN LINKBRIGHT™ (150 µl) is optimised for ~6 nmole of oligonucleotides.

*CPNs™ are available in a range of emission wavelengths (420 nm to 1130nm in the IR) see table or visit www.streambio.co.uk/products/

Procedure Overview

(Hands-on) – Reduce thiol oligonucleotides

Reaction time 30 min

(Hands-on) - Add oligonucleotides to a vial of LINKBRIGHT™

Reaction time 30 min

Stop reaction

(Hands-on) - Add Solution SC for 5 min

Before you Start

IMPORTANT: The purified oligonucleotide must be in a buffer free of ammonium ions, primary amines, or sodium azide preservatives, as they will disrupt the linkage reaction with the CPNs™. If the oligonucleotide is in, or has been lyophilized from an unsuitable buffer (e.g. Tris or glycine) or purified with ammonium sulphate, the buffer needs to be replaced with HEPES. Oligonucleotides can be purified and resuspended using standard methods, e.g. microdialysis or column separation.

Biomolecule Conjugation Procedure

Note: LINKBRIGHT™ Thiol Oligo Conjugation Kit is optimised for conjugation of oligonucleotides to the thiol group, for other biomolecules please either use an alternative kit, consult our guide, or contact technical support. LINKBRIGHT™ IgG antibody Conjugation Kits are also available

CPN Conjugation

1. Dissolve the lyophilized thiol-modified oligonucleotides in ddH₂O to the final concentration at 100 µM.
2. Take 60 µl of 100 µM thiol-modified oligonucleotide and add 60 µl of Solution HT (vial #2, pink cap)
3. Incubate at room temperature for 30 min
4. Add 150 µl of oligonucleotide[§] mixture to a vial of LINKBRIGHT™ CPN - Thiol
5. Incubate at room temperature for 30 min
6. Add 6µl Solution SC (vial #3, orange cap)
7. Incubate at room temperature for 5 min

(Optional) A magnetic separation[‡] method or a 100 kD MWCO desalting column can be chosen to purify the conjugates when required

[§]Optimal oligonucleotide to CPN ratio to be determined by end user.

[‡]CPNs™ can be attracted to magnets allowing the purifying and separating of the CPNs™ from unlinked reagents.